

DCUN1D1 antibody



Catalog Number: 109763

Product name

DCUN1D1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCUN1D1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, mouse kidney tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of DCUN1D1 recombinant protein (Accession Number: NM_020640). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109763(DCUN1D1 antibody) at dilution of



1:300

Immunohistochemical of paraffin-embedded human gliomas using Catalog No:109763(DCUN1D1 antibody) at dilution of 1:50 (under 10x lens)