

# DDHD2 antibody



Catalog Number: 109772

## General Information

### Product name:

DDHD2 antibody

### Immunogen:

[Human DDHD2 Recombinant protein \(GST tag\)](#)

### Applications:

ELISA, WB, IF

### Specificity:

Human, Mouse, Rat; other species not tested.

### Antibody description:

DDHD2 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in rat brain tissue, mouse lung tissue, mouse testis tissue. Observed molecular weight by Western-blot: 80-90 kDa

## Properties

### Preparation:

This antibody was obtained by immunization of DDHD2 recombinant protein (Accession Number: BC010504). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal

**Ig Type:**Rabbit IgG


## Applications

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109772(DDHD2 Antibody) at dilution of 1:2000



Immunofluorescent analysis of HEK-293 cells using Catalog No:109772(DDHD2 Antibody) at dilution of 1:25 and Alexa Fluor 594-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)