

# DCBLD2 antibody



Catalog Number: 109882

## Product name

DCBLD2 antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DCBLD2 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in human liver tissue, human brain tissue. Observed molecular weight by Western-blot: 93 kDa, 127 kDa

## Preparation

This antibody was obtained by immunization of DCBLD2 recombinant protein (Accession Number: BC029658). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109882(DCBLD2 antibody) at dilution of 1:500



IP Result of anti-DCBLD2 (IP:Catalog No:109882, 4ug; Detection:Catalog No:109882 1:500) with HepG2 cells lysate 2400ug.