

DCBLD2 antibody



Catalog Number: 109882

Product name

DCBLD2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DCBLD2 Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in human liver tissue, human brain tissue. Observed molecular weight by Western-blot: 93 kDa, 127 kDa

Preparation

This antibody was obtained by immunization of DCBLD2 recombinant protein (Accession Number: BC029658). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109882(DCBLD2 antibody) at dilution of 1:500



IP Result of anti-DCBLD2 (IP:Catalog No:109882, 4ug; Detection:Catalog No:109882 1:500) with HepG2 cells lysate 2400ug.