

# DENR antibody



Catalog Number: 109890

## Product name

DENR antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DENR Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells, A431 cells. Positive WB detected in A431 cells, HEK-293 cells, mouse brain tissue, mouse lung tissue. Observed molecular weight by Western-blot: 27-30kd

## Preparation

This antibody was obtained by immunization of DENR recombinant protein (Accession Number: NM\_003677). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109890(DENR antibody) at dilution of 1:1000



Immunofluorescent analysis of HEK-293 cells using Catalog No:109890(DENR Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)