

# DLEU1 antibody



Catalog Number: 109989

## Product name

DLEU1 antibody

## Specificity

human; other species not tested.

## Antibody description

DLEU1 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive IHC detected in human prostate cancer. Observed molecular weight by Western-blot: 27-30kd

## Preparation

This antibody was obtained by immunization of DLEU1 recombinant protein (Accession Number: BC020692). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109989(DLEU1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:109989(DLEU1 antibody) at dilution of 1:50 (under 40x lens)