

DRD1 antibody



Catalog Number: 110015

Product name

DRD1 antibody

Immunogen

[Human DRD1 Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

DRD1 Rabbit Polyclonal antibody. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of DRD1 recombinant protein (Accession Number: NM_000794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:110015(DRD1 antibody) at dilution of 1:1000