

# DRP1 antibody

Catalog Number: 110024

## Product name

DRP1 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DRP1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, NIH/3T3 cells, Transfected HEK-293 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 75-80 kDa

## Preparation

This antibody was obtained by immunization of DRP1 recombinant protein (Accession Number: NM\_012063). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110024(DNM1L,DLP1 Antibody) at dilution of



IP Result of anti-DNM1L, DLP1 (IP:Catalog No:110024, 4ug; Detection:Catalog No:110024 1:500) with mouse brain tissue lysate 3600ug.