

DTWD2 antibody

Catalog Number: 110035

Product name

DTWD2 antibody

Specificity

Human; other species not tested.

Antibody description

DTWD2 Rabbit Polyclonal antibody. Positive IP detected in COLO 320 cells. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of DTWD2 recombinant protein (Accession Number: NM_173666). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IP

Dilutions

Recommended Dilution:

IP: 1:200-1:1000

Validations



IP Result of anti-DTWD2 (IP:Catalog No:110035, 3ug; Detection:Catalog No:110035 1:300) with COLO 320 cells lysate 2000ug.