

DTX1 antibody

Catalog Number: 110036

Product name

DTX1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DTX1 Rabbit Polyclonal antibody. Positive IF detected in mouse olfactory epithelium tissue. Positive IHC detected in human testis tissue.

Preparation

This antibody was obtained by immunization of DTX1 recombinant protein (Accession Number: XM_011538009). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, IHC


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: N/A

Validations



Fluorescent IHC on 1%PLP fixed frozen mouse olfactory epithelium tissue of DTX1 antibody (Catalog No:110036). Basal cell staining as expected with Sus cell staining on the apical layer. By Dr. Brian Lin (Schwob Lab).



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110036(DTX1 Antibody) at dilution of 1:50 (under 40x lens)