

DOCK6 antibody



Catalog Number: 110045

Product name

DOCK6 antibody

Immunogen

[Human DOCK6 Recombinant protein \(His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

DOCK6 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells. Observed molecular weight by Western-blot: 230 kDa

Preparation

This antibody was obtained by immunization of DOCK6 recombinant protein (Accession Number: NM_020812). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110045(Dock6 Antibody) at dilution of 1:300