

# DOM3Z antibody

Catalog Number: 110056

## Product name

DOM3Z antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DOM3Z Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, SH-SY5Y cells. Positive IP detected in mouse testis tissue. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of DOM3Z recombinant protein (Accession Number: XM\_047418262). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000


## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110056(DOM3Z antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110056(DOM3Z antibody) at dilution of 1:50



IP Result of anti-DOM3Z (IP:Catalog No:110056,

# DOM3Z antibody



Catalog Number: 110056

---

4ug; Detection:Catalog No:110056 1:300) with mouse testis tissue lysate 4000ug.