

# DNAJB3 antibody

Catalog Number: 110094

## Product name

DNAJB3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DNAJB3 Rabbit Polyclonal antibody. Positive IHC detected in human testis. Positive WB detected in mouse testis tissue, Jurkat cells, PC-3 cells. Observed molecular weight by Western-blot: 28-30kd

## Preparation

This antibody was obtained by immunization of DNAJB3 recombinant protein (Accession Number: NM\_001001394). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110094(DNAJB3 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human testis using Catalog No:110094(DNAJB3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:110094(DNAJB3 antibody) at dilution of 1:100 (under 40x lens)