

DNAJB7 antibody



Catalog Number: 110098

Product name

DNAJB7 antibody

Immunogen

[Human DNAJB7 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DNAJB7 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of DNAJB7 recombinant protein (Accession Number: NM_145174). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110098(DNAJB7 Antibody) at dilution of 1:300