

DNAJC10 antibody



Catalog Number: 110102

Product name

DNAJC10 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

DNAJC10 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse liver tissue. Positive IHC detected in human pancreas cancer. Observed molecular weight by Western-blot: 75-80kd

Preparation

This antibody was obtained by immunization of DNAJC10 recombinant protein (Accession Number: BC034713). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110102(DNAJC10 antibody) at dilution of 1:1500



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110102(DNAJC10 antibody) at dilution of 1:100 (under 10x lens)