

# DNAJC15 antibody



Catalog Number: 110105

## Product name

DNAJC15 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DNAJC15 Rabbit Polyclonal antibody. Positive IHC detected in human liver. Positive WB detected in mouse liver tissue, human brain tissue. Observed molecular weight by Western-blot: 18-20kd

## Preparation

This antibody was obtained by immunization of DNAJC15 recombinant protein (Accession Number: NM\_013238). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110105(DNAJC15 antibody) at dilution of



1:1000  
Immunohistochemical of paraffin-embedded human liver using Catalog No:110105(DNAJC15 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:110105(DNAJC15 antibody) at dilution of 1:100 (under 40x lens)