

# DNAJC16 antibody



Catalog Number: 110106

## Product name

DNAJC16 antibody

## Immunogen

[Human DNAJC16 Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DNAJC16 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 91kd

## Preparation

This antibody was obtained by immunization of DNAJC16 recombinant protein (Accession Number: NM\_015291). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110106(DNAJC16 antibody) at dilution of 1:1500