

# DUS4L antibody



Catalog Number: 110113

## Product name

DUS4L antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DUS4L Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue. Observed molecular weight by Western-blot: 22,30 kDa

## Preparation

This antibody was obtained by immunization of DUS4L recombinant protein (Accession Number: BC111774). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110113(DUS4L antibody) at dilution of 1:1000