

DYNC2LI1 antibody



Catalog Number: 110139

Product name

DYNC2LI1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DYNC2LI1 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, human kidney tissue. Observed molecular weight by Western-blot: 40 kDa

Preparation

This antibody was obtained by immunization of DYNC2LI1 recombinant protein (Accession Number: NM_001193464). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110139(DYNC2LI1 antibody) at dilution of 1:500