

# DYRK1B antibody



Catalog Number: 110144

## Product name

DYRK1B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DYRK1B Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HeLa cells, human colon tissue, human heart tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 57 kDa

## Preparation

This antibody was obtained by immunization of DYRK1B recombinant protein (Accession Number: XM\_005259398). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110144(DYRK1B antibody) at dilution of 1:500