

Dysbindin antibody



Catalog Number: 110147

Product name

Dysbindin antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Dysbindin Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human testis tissue, mouse testis tissue. Observed molecular weight by Western-blot: 48kd, 36kd

Preparation

This antibody was obtained by immunization of Dysbindin recombinant protein (Accession Number: NM_032122). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110147(DTNBP1 antibody) at dilution of 1:1000