

DYX1C1 antibody

Catalog Number: 110150

Product name

DYX1C1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

DYX1C1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, rat testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of DYX1C1 recombinant protein (Accession Number: NM_001033560). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110150(DYX1C1 antibody) at dilution of



1:1000

IP Result of anti-DYX1C1 (IP:Catalog No:110150, 4ug; Detection:Catalog No:110150 1:500) with mouse testis tissue lysate 4000ug.