

# EHHADH antibody



Catalog Number: 110227

## Product name

EHHADH antibody

## Immunogen

[Human EHHADH Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

EHHADH Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue, SMMC-7721 cells. Observed molecular weight by Western-blot: 69 kDa

## Preparation

This antibody was obtained by immunization of EHHADH recombinant protein (Accession Number: NM\_001966). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:1000-1:10000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110227(EHHADH Antibody) at dilution of 1:2000