

# EHMT2/G9a antibody



Catalog Number: 110228

## Product name

EHMT2/G9a antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EHMT2/G9a Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 140 kDa

## Preparation

This antibody was obtained by immunization of EHMT2/G9a recombinant protein (Accession Number: NM\_001318833). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110228(EHMT2 antibody) at dilution of 1:400