

EHMT2/G9a antibody



Catalog Number: 110228

General Information

Product name:

EHMT2/G9a antibody

Immunogen:

Applications:

ELISA, WB

Specificity:

Human, Mouse, Rat; other species not tested.

Antibody description:

EHMT2/G9a Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 140 kDa

Properties

Preparation:

This antibody was obtained by immunization of EHMT2/G9a recombinant protein (Accession Number: NM_001318833). Purification method: Antigen affinity purified.

Formulation:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage:

Store at -20°C. DO NOT ALIQUOT

Clonality:

Polyclonal


Ig Type: Rabbit IgG

Applications

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110228(EHMT2 antibody) at dilution of 1:400