

EMILIN2 antibody



Catalog Number: 110240

Product name

EMILIN2 antibody

Specificity

Human; other species not tested.

Antibody description

EMILIN2 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue.

Preparation

This antibody was obtained by immunization of EMILIN2 recombinant protein (Accession Number: NM_032048). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:110240(EMILIN2 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using Catalog No:110240(EMILIN2 Antibody) at dilution of 1:200 (under 40x lens).