

# EDG2 antibody

Catalog Number: 110303

## Product name

EDG2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EDG2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, A375 cells, HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 41-50 kDa

## Preparation

This antibody was obtained by immunization of EDG2 recombinant protein (Accession Number: NM\_001387497). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



WB result of anti-LPAR1 (Catalog No.:1:1,000).



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110303(LPAR1,EDG2 antibody) at dilution of 1:1000