

EPB41 antibody

Catalog Number: 110343

Product name

EPB41 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EPB41 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, Jurkat cells. Positive IP detected in HL-60 cells. Observed molecular weight by Western-blot: 80kd, 135kd

Preparation

This antibody was obtained by immunization of EPB41 recombinant protein (Accession Number: NM_001166006). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:110343(EPB41 antibody) at dilution of 1:100



IP Result of anti-EPB41 (IP:Catalog No:110343, 4ug; Detection:Catalog No:110343 1:200) with HL-60 cells lysate 4000ug.