

# EXOSC8 antibody



Catalog Number: 110414

## Product name

EXOSC8 antibody

## Immunogen

[Human EXOSC8 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

EXOSC8 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human colon tissue, human pancreas cancer tissue. Positive IF detected in 4T1 cells. Positive WB detected in HeLa cells, HEK-293 cells, human brain tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 30-35 kDa

## Preparation

This antibody was obtained by immunization of EXOSC8 recombinant protein (Accession Number: NM\_181503). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF, IHC, IP

## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000


IHC: 1:20-1:200

IF: 1:20-1:200

## Validations




HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110414(EXOSC8 antibody) at dilution of 1:500




Immunofluorescent analysis of 4T1 cells, using EXOSC8 antibody Catalog No:110414 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent

# EXOSC8 antibody


Catalog Number: 110414



IP Result of anti-EXOSC8 (IP:Catalog No:110414, 4ug; Detection:Catalog No:110414 1:500) with HEK-293 cells lysate 1560ug.



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110414 (EXOSC8 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:110414 (EXOSC8 Antibody) at dilution of 1:50 (under 40x lens)