

F8 antibody



Catalog Number: 110432

Product name

F8 antibody

Specificity

Human; other species not tested.

Antibody description

F8 Rabbit Polyclonal antibody. Positive IF detected in human umbilical vein endothelial cells, HUVEC cells. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of F8 recombinant protein (Accession Number: NM_019863). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: N/A

Validations



IF result of anti-F8(Catalog No:110432) with primary human umbilical vein endothelial cells.



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110432(F8 Antibody) at dilution of 1:300