

# FABP7 antibody



Catalog Number: 110442

## General Information

### Product name:

FABP7 antibody

### Immunogen:

### Applications:

ELISA, WB

### Specificity:

Human, Mouse, Rat; other species not tested.

### Antibody description:

FABP7 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 15kd

## Properties

### Preparation:

This antibody was obtained by immunization of FABP7 recombinant protein (Accession Number: BC012299). Purification method: Antigen affinity purified.

### Formulation:

PBS with 0.02% sodium azide and 50% glycerol

pH 7.3.

### Storage:

Store at -20°C. DO NOT ALIQUOT

### Clonality:

Polyclonal


**Ig Type:** Rabbit IgG

## Applications

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110442(FABP7 antibody) at dilution of 1:800