

FADS3 antibody



Catalog Number: 110450

Product name

FADS3 antibody

Immunogen

[Human FADS3 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

FADS3 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, human brain tissue, mouse kidney tissue. Observed molecular weight by Western-blot: 49-51kd

Preparation

This antibody was obtained by immunization of FADS3 recombinant protein (Accession Number: NM_021727). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:110450(FADS3 antibody) at dilution of 1:300