

# FAH antibody

Catalog Number: 110453

---

## Product name

FAH antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FAH Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 40-45kd

## Preparation

This antibody was obtained by immunization of FAH recombinant protein (Accession Number: NM\_000137). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110453(FAH antibody) at dilution of 1:500