

FAM134B antibody



Catalog Number: 110473

Product name

FAM134B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FAM134B Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, mouse brain tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of FAM134B recombinant protein (Accession Number: BC030517). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110473(FAM134B Antibody) at dilution of 1:500