

ESR2 Beta-3 antibody

Catalog Number: 110475

Product name

ESR2 Beta-3 antibody

Specificity

Human; other species not tested.

Antibody description

ESR2 Beta-3 Rabbit Polyclonal antibody. Positive IHC detected in human cervix tissue, human breast cancer tissue, human endometrial cancer tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001437). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human cervix using Catalog No:110475(ESR2 Beta-3 antibody) at dilution of 1:50 (under 10x lens)