

EXOC6 antibody

Catalog Number: 110503

Product name

EXOC6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

EXOC6 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 116 kDa

Preparation

This antibody was obtained by immunization of EX^oC6 recombinant protein (Accession Number: NM_019053). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110503(EXoc6 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using EXoc6 antibody Catalog No:110503 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).