

FAM98C antibody



Catalog Number: 110513

Product name

FAM98C antibody

Specificity

Human; other species not tested.

Antibody description

FAM98C Rabbit Polyclonal antibody. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 37 kDa; 28 kDa

Preparation

This antibody was obtained by immunization of FAM98C recombinant protein (Accession Number: BC034383). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:110513(FAM98C Antibody) at dilution of 1:300