

# FAU antibody

Catalog Number: 110534

## Product name

FAU antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FAU Rabbit Polyclonal antibody. Positive IHC detected in human stomach tissue, human colon cancer tissue. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 14kd

## Preparation

This antibody was obtained by immunization of FAU recombinant protein (Accession Number: NM\_001997). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110534(FAU antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:110534(FAU Antibody) at dilution of 1:50 (under 40x lens)