

# FBXO43 antibody



Catalog Number: 110555

## Product name

FBXO43 antibody

## Specificity

Human; other species not tested.

## Antibody description

FBXO43 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Positive IF detected in HEK-293 cells. Observed molecular weight by Western-blot: ~70 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001029860). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110555(FBXO43 Antibody) at dilution of 1:300



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using Catalog No:110555(FBXO43 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)