

# FCHSD1 antibody



Catalog Number: 110568

---

## Product name

FCHSD1 antibody

## Specificity

human; other species not tested.

## Antibody description

FCHSD1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 46-50 kDa

## Preparation

This antibody was obtained by immunization of FCHSD1 recombinant protein (Accession Number: BC047016). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110568(FCHSD1 antibody) at dilution of 1:1000