

# FAM186B antibody



Catalog Number: 110581

## Product name

FAM186B antibody

## Specificity

Human; other species not tested.

## Antibody description

FAM186B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 120 kd

## Preparation

This antibody was obtained by immunization of FAM186B recombinant protein (Accession Number: XM\_006719625). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110581(FAM186B Antibody) at dilution of 1:600