

FAM40B antibody



Catalog Number: 110588

Product name

FAM40B antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FAM40B Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, HeLa cells, HepG2 cells, Jurkat cells, K-562 cells, MCF-7 cells, mouse colon tissue, mouse testis tissue, Raji cells. Observed molecular weight by Western-blot: 95 kDa

Preparation

This antibody was obtained by immunization of FAM40B recombinant protein (Accession Number: NM_001134336). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:110588(FAM40B Antibody) at dilution of 1:600