

# FAM62A antibody



Catalog Number: 110595

---

## Product name

FAM62A antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

FAM62A Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells, mouse uterus tissue. Observed molecular weight by Western-blot: 125-130kd,140-150kd

## Preparation

This antibody was obtained by immunization of FAM62A recombinant protein (Accession Number: NM\_015292). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110595(FAM62A antibody) at dilution of 1:500