

# FDPS antibody



Catalog Number: 110608

## Product name

FDPS antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FDPS Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HepG2 cells, human cerebellum tissue, human spleen tissue, MCF7 cells, mouse liver tissue, rat liver tissue. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 39-41 kDa

## Preparation

This antibody was obtained by immunization of FDPS recombinant protein (Accession Number: NM\_002004). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110608(FDPS antibody) at dilution of 1:1500



IP Result of anti-FDPS (IP:Catalog No:110608, 4ug; Detection:Catalog No:110608 1:1000) with mouse liver tissue lysate 4000ug.