

# FKBP1A antibody

Catalog Number: 110671

## Product name

FKBP1A antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FKBP1A Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human heart tissue. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 12-16 kDa

## Preparation

This antibody was obtained by immunization of FKBP1A recombinant protein (Accession Number: NM\_054014). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:110671(FKBP1A antibody) at dilution of 1:100



(under 25x lens)

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110671(FKBP1A Antibody) at dilution of 1:1000