

# FKBPL antibody

Catalog Number: 110684

## Product name

FKBPL antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FKBPL Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human breast cancer tissue. Positive IP detected in MCF-7 cells. Positive WB detected in MCF7 cells, COLO 320 cells, HEK-293 cells, HeLa cells, human brain tissue, human testis tissue, Jurkat cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 42-48 kDa; 38-40 kDa

## Preparation

This antibody was obtained by immunization of FKBPL recombinant protein (Accession Number: NM\_022110). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110684(FKBPL antibody) at dilution of 1:500




IP Result of anti-FKBPL (IP:Catalog No:110684, 3ug; Detection:Catalog No:110684 1:800) with MCF-7 cells lysate 2500ug.


# FKBPL antibody

Catalog Number: 110684

---



Immunohistochemical of paraffin-embedded human testis using Catalog No:110684(FKBPL antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:110684(FKBPL antibody) at dilution of 1:50 (under 40x lens)