

Flag tag antibody

Catalog Number: 110685

Product name

Flag tag antibody

Specificity

Tag

Antibody description

Flag tag Rabbit Polyclonal antibody. Positive IF detected in Transfected HEK-293 cells. Positive WB detected in recombinant protein.

Preparation

This antibody was obtained by immunization of Flag tag recombinant protein. Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:4000

IF: 1:10-1:100

Validations



Western blot of Flagged fusion protein with anti-Flag (Catalog No:117313) at various dilutions.



Immunofluorescent analysis of Transfected HEK-293 cells using Catalog No:117313(Flag tag Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)