

FLI1 antibody

Catalog Number: 110687

Product name

FLI1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FLI1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in HL-60 cells, human brain tissue, Jurkat cells, mouse thymus tissue. Observed molecular weight by Western-blot: 51 kDa

Preparation

This antibody was obtained by immunization of FLI1 recombinant protein (Accession Number: NM_002017). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:50-1:500


Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:110687(FLI1 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:110687(FLI1 Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog

FLI1 antibody



Catalog Number: 110687

No:110687(FLI1 Antibody) at dilution of 1:200

(under 40x lens). heat mediated antigen retrieved with Tris-EDTA buffer(pH9).