

FOXC1 antibody



Catalog Number: 110692

Product name

FOXC1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FOXC1 Rabbit Polyclonal antibody. Positive WB detected in MCF-7 cells, HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001453). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110692(FOXC1 Antibody) at dilution of 1:300