

FOXR1 antibody

Catalog Number: 110714

Product name

FOXR1 antibody

Immunogen

[Human FOXR1 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

FOXR1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in fetal human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of FOXR1 recombinant protein (Accession Number: NM_181721). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110714(FOXR1 Antibody) at dilution of 1:600



IP Result of anti-FOXR1 (IP:Catalog No:110714, 4ug; Detection:Catalog No:110714 1:500) with mouse brain tissue lysate 4000ug.